




Everylab sample preparation packages

Based on your laboratory needs, benefit from our Everylab packages with either a ventilated or refrigerated microcentrifuge and a manual pipette kit, as well as tubes.



Thermo Scientific™ Pico™ and Fresco™ 21 Microcentrifuges

- **Ease of use**—Simplify operation with intuitive controls, easy-to-read displays and fast one-click centrifuge lid closure
- **Efficient**—Lightweight rotors are designed to provide acceleration up to 21,100 x g in 12 seconds. Outstanding corrosion resistance with polymer rotors; includes a broad array of microcentrifuge rotors and accessories for processing micro-volumes
- **Ergonomic**—One-click biocontainment lid opening and closing, as well as transparent design for added safety with certified¹ Thermo Scientific™ ClickSeal™ Biocontainment Lids



Transparent ClickSeal biocontainment lids



Fast one-click centrifuge lid closure



Large, bright and easy-to-read display



Up to 21,100 x g



Thermo Scientific™ FinnpiPETTE™ F2 Pipettes

- **Rugged and durable**—The FinnpiPETTE F2 system contains tough polyvinylidene fluoride (PVDF) components that stand up to harsh chemicals and the damaging effects of UV light; its rugged design withstands physical use without damage
- **Fully autoclavable**— Decontamination of the FinnpiPETTE F2 is simple because there's no need to disassemble and reassemble the pipette for autoclaving, helping minimize disruption and downtime
- **Maintenance is easy**—Simply detach the tip cone for efficient daily maintenance or decontamination when not using an autoclave

Ordering information

Description	Cat. No.
Everylab package with refrigerated centrifuge*	76906-206
Thermo Scientific™ Fresco™ 21 Centrifuge, 230 V incl. 24 x 1.5/2 mL rotor with ClickSeal biocontainment lid Thermo Scientific™ FinnpiPETTE™ F2 GLP Kit 2: <ul style="list-style-type: none"> • Single-channel 0.2-2 µL, 2-20 µL, 20-200 µL and 100-1,000 µL pipettes • Flex tips: 10 micro: 1 x 96; 200: 2 x 96; 1000: 1 x 96 • F-stand • Good Laboratory Pipetting Guide • Reagent reservoir demo pack and • 1.5 Snap Cap Low Retention microcentrifuge tubes, full case (5,000 tubes) 	
Everylab package with ventilated centrifuge*	76906-204
Thermo Scientific™ Pico™ 21 Centrifuge, 230 V incl. 24 x 1.5/2 mL rotor with ClickSeal biocontainment lid Thermo Scientific™ FinnpiPETTE™ F2 GLP Kit 2: <ul style="list-style-type: none"> • Single-channel 0.2-2 µL, 2-20 µL, 20-200 µL and 100-1,000 µL pipettes • Flex tips: 10 micro: 1 x 96; 200: 2 x 96; 1000: 1 x 96 • F-stand • Good Laboratory Pipetting Guide • Reagent reservoir demo pack and • 1.5 Snap Cap Low Retention microcentrifuge tubes, full case (5,000 tubes) 	

1



2



*This package is not considered a medical device and has not been evaluated for use in environments or applications involving the diagnosis of diseases or other conditions, or in the cure, mitigation, treatment, or prevention of disease, in man or other animals.