



Technical Data Sheet




VWR® Cell Lifter, 19mm Flat Blade with 9mm J-Hook end

Description: Cell Lifter, 19mm Flat Blade with 9mm J-Hook end, Sterilized ,Non-pyrogenic


Purpose: Flat Blade Cell Lifter is manufactured with high-grade polyethylene(PE) which is tensile to ensure minimal cell damage during manual harvesting

- Materials
- 19mm Flat Blade: HDPE/High-density Polyethylene    Color: Green
  - 9mm J-Hook end: HDPE/High-density Polyethylene    Color: Green



- Features
- Convenient operation,special blade design minimizes cell damage
  - Individually wrapped
  - DNase/RNase free and Non-pyrogenic
  - Sterilized by irradiation SAL10<sup>-6</sup> (ISO11137)
  - Shelf life: 3 years after month of production
  - Manufactured in a class 100,000 room environment
  - Manufactured under ISO13485 and ISO9001 quality management system

Dimensions ( Unit: mm)



VWR® Cell Lifter, 19mm Flat Blade with 9mm J-Hook end

VWR NA ENA. No.	VWR EU Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
76036-004	734-2978	19mm Flat Blade with 9mm J-Hook end, Sterilized ,Non-pyrogenic	Yes	1	100	100



Technical Data Sheet



VWR® Cell Lifter, 19mm Flat Blade with 2.5mm Narrow Blade end

Description: Cell Lifter, 19mm Flat Blade with 9mm J-Hook end, Sterilized ,Non-pyrogenic

Purpose: Flat Blade Cell Lifter is manufactured with high-grade polyethylene(PE) which is tensile to ensure minimal cell damage during manual harvesting


Materials

19mm Flat Blade: HDPE/High-density Polyethylene    Color: Green


2.5mm Narrow Blade end: HDPE/High-density Polyethylene    Color: Green

Features

- Convenient operation, special blade design minimizes cell damage
- Individually wrapped
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL10<sup>-6</sup> (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



Dimensions ( Unit: mm)



VWR® Cell Lifter, 19mm Flat Blade with 2.5mm Narrow Blade end

VWR NA ENA. No.	VWR EU Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
76036-006	734-2979	19mm Flat Blade with 2.5mm Narrow Blade end, Sterilized ,Non-pyrogenic	Yes	1	100	100