



## Interior

Front View



Side View



Cat# 76284-668/MFS-15-C



| Cu.Ft. | Defrost | Door    | Int Door | Shelves   | Interior Dimensions |        |        | Refrigerant | H.P. | Amps | Weight   |
|--------|---------|---------|----------|-----------|---------------------|--------|--------|-------------|------|------|----------|
|        |         |         |          |           | W"                  | D"     | H"     |             |      |      |          |
| 15     | Manual  | 1 Solid | 0        | 2 Baskets | 49 3/8              | 21 1/2 | 25 5/8 | R600a       | 1/4  | 5    | 137 lbs. |

Note: Interior dimensions are estimates and subject to change without notice