





External Threaded Self Standing / Star foot / Round Bottom Vials



Internal Threaded Self Standing / Star foot / Round Bottom vials

VWR Global Cat#	Description	Capacity ml	Brim Volume ml	Thread type	Foot type	Neck Internal dia	Tube OD (mm)	General Wall thickness (mm)	Height with Closure (mm)	Closure OD (mm)	Weight (gms)	Graduation range	Gamma Sterilized
VWRI210-000001	CRYOVIAL EXT TRD SF 1/1.2 ML	1 / 1.2ml	1.8ml	External Threaded	Star Foot	9.8	12	1.2	35	13.8	1.53	0.5 - 1ml	Yes
VWRI210-000003	CRYOVIAL EXT TRD SF 1.8/2 ML	1.8 / 2ml	2.6ml	External Threaded	Star Foot	9.8	12	1.2	46.9	13.8	1.98	0.5 - 1.8ml	Yes
VWRI210-000007	CRYOVIAL EXT TRD SF 3 ML	3ml	4.2ml	External Threaded	Star Foot	9.8	12	1.2	69.6	13.8	2.79	0.5 - 3.0ml	Yes
VWRI210-000008	CRYOVIAL EXT TRD SF 3.6/4 ML	3.6 / 4ml	4.5ml	External Threaded	Star Foot	9.8	12.1	1.2	74.1	13.8	2.95	0.5 - 3.0ml	Yes
VWRI210-000011	CRYOVIAL EXT TRD SF 4.5/5 ML	4.5 / 5ml	5.5ml	External Threaded	Star Foot	9.8	12	1.2	89.9	13.8	3.38	0.5 - 4.5ml	Yes
VWRI210-000002	CRYOVIAL EXT TRD SS 1/1.2 ML	1 / 1.2ml	1.8ml	External Threaded	Self Standing	9.8	12.1	1.2	41.2	13.8	1.91	0.5 - 1ml	Yes
VWRI210-000004	CRYOVIAL EXT TRD SS 1.8/2 ML	1.8 / 2ml	2.6ml	External Threaded	Self Standing	9.8	12	1.22	48	13.8	2.02	0.5 - 1.8ml	Yes
VWRI210-000006	CRYOVIAL EXT TRD SS 3 ML	3ml	4.2ml	External Threaded	Self Standing	9.8	12.1	1.2	70.7	13.8	2.86	0.5 - 3ml	Yes
VWRI210-000009	CRYOVIAL EXT TRD SS 3.6/4 ML	3.6 / 4ml	4.5ml	External Threaded	Self Standing	9.8	12.1	1.2	75.2	13.8	3.02	0.5 - 3.5ml	Yes
VWRI210-000010	CRYOVIAL EXT TRD SS 4.5/5 ML	4.5 / 5ml	5.5ml	External Threaded	Self Standing	9.8	12	1.22	91	13.8	3.44	0.5 - 4.5ml	Yes
VWRI210-000005	CRYOVIAL EXT TRD RB 1.8/2 ML	1.8 / 2ml	2.6ml	External Threaded	Round Bottom	9.8	12.1	1.2	46	13.8	2	0.5 - 1.8ml	Yes
VWRI210-000012	CRYOVIAL INT TRD SF 1/1.2 ML W/WASHER	1 / 1.2ml	1.2ml	Internal Threaded	Star Foot	10.4	12.4	1.2	41.5	12.2	1.81	0.5 - 1ml	Yes
VWRI210-000013	CRYOVIAL INT TRD SF 1.8/2 ML W/WASHER	1.8 / 2ml	2.2ml	Internal Threaded	Star Foot	10.4	12.4	1.05	47.8	12.2	1.88	0.5 - 1.8ml	Yes
VWRI210-000017	CRYOVIAL INT TRD SF 3.6/4 ML W/WASHER	3.6 / 4ml	5ml	Internal Threaded	Star Foot	10.4	12.4	1	77.5	12.2	2.8	0.5 - 3.5ml	Yes
VWRI210-000019	CRYOVIAL INT TRD SF 4.5/5 ML W/WASHER	4.5 / 5ml	5.5ml	Internal Threaded	Star Foot	10.4	12.4	1.15	90.3	12.2	3.4	0.5 - 4.5ml	Yes
VWRI210-000014	CRYOVIAL INT TRD SS 1.8/2 ML W/WASHER	1.8 / 2ml	2.2ml	Internal Threaded	Self Standing	10.4	12.4	1.05	48.9	12.2	1.93	0.5 - 1.8ml	Yes
VWRI210-000016	CRYOVIAL INT TRD SS 3.6/4 ML W/WASHER	3.6 / 4ml	5.0ml	Internal Threaded	Self Standing	10.4	12.4	1	78.6	12.2	2.89	0.5 - 3.5ml	Yes
VWRI210-000020	CRYOVIAL INT TRD SS 4.5/5 ML W/WASHER	4.5 / 5ml	6.2ml	Internal Threaded	Self Standing	10.4	12.4	1	92.8	12.2	2.69	0.5 - 4.5ml	Yes
VWRI210-000015	CRYOVIAL INT TRD RB 1.8/2 ML W/WASHER	1.8 / 2ml	2.8ml	Internal Threaded	Round Bottom	10.4	12.4	1	49.2	12.2	1.9	0.5 - 1.8ml	Yes
VWRI210-000018	CRYOVIAL INT TRD RB 3.6/4 ML W/WASHER	3.6 / 4ml	5.0ml	Internal Threaded	Round Bottom	10.4	12.4	1	77	12.2	2.78	0.5 - 2ml	Yes
VWRI210-000021	CRYOVIAL INT TRD RB 4.5/5 ML W/WASHER	4.5 / 5ml	5.5ml	Internal Threaded	Round Bottom	10.4	12.4	1	91	12.2	3.22	0.5 - 4.5ml	Yes
VWRI210-000026	CRYOVIAL 2D EXT TRD W/CAP SS 2ML	2ml	2.6ml	External Threaded	Self Standing	9.8	12	1.22	48	13.8	2.02	0.5 - 1.8ml	Yes
VWRI210-000027	CRYOVIAL 2D IN 81 RACK BLUE 2ML	2ml	2.6ml	External Threaded	Self Standing	9.8	12	1.22	48	13.8	2.02	0.5 - 1.8ml	Yes
VWRI210-000028	CRYOVIAL 2D IN 81 RACK GREEN 2ML	2ml	2.6ml	External Threaded	Self Standing	9.8	12	1.22	48	13.8	2.02	0.5 - 1.8ml	Yes
VWRI210-000029	CRYOVIAL 2D IN 81 RACK RED 2ML	2ml	2.6ml	External Threaded	Self Standing	9.8	12	1.22	48	13.8	2.02	0.5 - 1.8ml	Yes
VWRI210-000030	CRYOVIAL 2D IN 81 RACK YELLOW 2ML	2ml	2.6ml	External Threaded	Self Standing	9.8	12	1.22	48	13.8	2.02	0.5 - 1.8ml	Yes

Material

Vial : PP 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

Closure : HDPE, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Clean Statement: Free of detectable DNase, RNase & Pyrogen

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Application Temperature: -196°C to +121°C

Leak Proof :

Material Recycle Code



Description: Cryovials - Sterile

Drawing:
Not to Scale



Date: 03/28/2023

Revision: 0

Cryo Vial Product type

External Thread Vial



Internal Thread Vial



Foot Type

Self Standing(SS)



Star Foot (SF)



Round Bottom (RB)



Self Standing (SS) 2D Coded



Cryo Coder / Cap Insert



VWR Global Cat#	Description
VWRI210-000031	CAP INSERT VIOLET
VWRI210-000032	CAP INSERT WHITE
VWRI210-000033	CAP INSERT 5 COLORS
VWRI210-000034	CAP INSERT CRYOVIAL BLUE
VWRI210-000035	CAP INSERT CRYOVIAL BURTORANGE
VWRI210-000036	CAP INSERT CRYOVIAL LILAC
VWRI210-000037	CAP INSERT CRYOVIAL RED
VWRI210-000038	CAP INSERT CRYOVIAL TAN
VWRI210-000039	CAP INSERT CRYOVIAL YELLOW
VWRI210-000040	CAP INSERT GREEN
VWRI210-000041	CAP INSERT GREY
VWRI210-000042	CAP INSERTCRYOVIAL PINK

Description: Cryovials - Sterile

Drawing:
Not to Scale



Date: 03/28/2023

Revision: 0