

KIMTECH™

BUILT FOR
SCIENCE.
—
MADE
FOR YOU.

avantor™
delivered by **vwr™**

TEXTURED FINGERTIPS
TO ENHANCE WET AND
DRY GRIP

0.65 AQL
FOR PINHOLES

PROVIDES
PROTECTION AGAINST
HARSH CHEMICALS*

Kimtech™ Polaris™ Nitrile Exam Gloves

9.5"
length

12"
length

With protection, durability and
comfort you can count on,
these gloves are built for
science, made for you.

*Tested against 29 common chemicals,
Chemotherapy Drugs, Opioid Fentanyl
Citrate, Simulated Gastric Acid and
Fentanyl in Simulated Gastric Acid:
Per ASTM D 6978-05





INTRODUCING
Kimtech™ Polaris™
Nitrile Exam Gloves







COMFORT + PROTECTION



Kimtech™ is your partner in driving progress. As part of the 150-year-old Kimberly-Clark legacy of innovation and partnership with the scientific community, we're proud to present our next step into the future – Kimtech™ Polaris™ Nitrile Exam Gloves.

The proprietary high-quality nitrile formulation delivers the highest level of protection in the Kimtech™ exam glove portfolio. With protection, durability and comfort you can count on, these gloves are built for science, made for you.

Kimtech™ Polaris™ Nitrile Exam Gloves meet and exceed the industry standards for exam gloves.

-  Designed for use in laboratories; also suitable for pharmaceutical manufacturing, med device manufacturing, biotechnology.
-  Tested against 24 chemotherapy drugs, which is almost triple the number required by the ASTM D6978 standard for chemotherapy use.
-  More than double the strength of the industry standard, the Kimtech™ Polaris™ Nitrile Exam Glove has an average tensile strength of 31 Mpa (after aging).
-  Tested for use with Fentanyl and Simulated Gastric Acid per ASTM D6978-05.
-  Meets FDA food contact requirements set forth in 21 CFR 177.2600.
-  Does not contain silicone according to results obtained using test method SAE AMS 3819C.

Apex Make your laboratory exceptional.

The Apex™ Program is a complimentary training program that ensures a seamless transition to the Kimtech™ Brand of products in environments, such as research labs, semi-conductor production, or advanced manufacturing, where personal safety and contamination control are essential. By effectively engaging employees, Apex ensures the protection of your process while enhancing efficiency and performance.

*When used in non-hazardous applications as classified by local regulations.

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you for your initial purchase, via FREE product, for up to \$1,000.

™/® Trademarks of Kimberly-Clark Worldwide, Inc., or their respective holders and may be used under license. Marques de commerce de Kimberly-Clark Worldwide, Inc. ou de leurs détenteurs respectifs et peuvent être utilisées sous licence. © KCWW T2332K-VWR 1/24

Kimtech™ Polaris™ Nitrile Exam Gloves | 9.5"

SIZE	KCP CODE	VWR CODE	BOX QTY	CASE QTY
XS	62770	77308-536	100	1000
S	62771	77308-538	100	1000
M	62772	77308-540	100	1000
L	62773	77308-542	100	1000
XL	62774	77308-544	90	900

Kimtech™ Polaris™ Xtra Nitrile Exam Gloves | 12"

SIZE	KCP CODE	VWR CODE	BOX QTY	CASE QTY
XS	62760	77308-526	50	500
S	62761	77308-528	50	500
M	62762	77308-530	50	500
L	62763	77308-532	50	500
XL	62764	77308-534	50	500



These gloves are recyclable through our RightCycle™ Program

The RightCycle™ Program is a sustainability initiative by Kimberly-Clark Professional that is committed to leading the way to zero waste in landfills. The Kimtech™ Polaris™ Nitrile Exam Gloves are recyclable through our RightCycle™ Program, along with other PPE, such as protective clothing and safety goggles.*