



TRIO.BASTM

BIOLOGICAL AIR SAMPLER

Now With IQ, OQ, PQ Documents

Microbial Air Sampler Kits



order on
VWR.COM

MINI



Each kit comes with the following:

TRIO.BAS MINI BLUETOOTH air sampler, battery charger and cable, aspirating head, cover head, LIGHT carrying case, and calibration certificate.

*IQ, OQ, PQ documents not available with TRIO.BAS MINI kit

VWR Cat. no.

MINI Kit, 100 liters/min., contact plate	76076-694
MINI Kit, 100 liters/min., Petri plate	76169-638
MINI Kit, 200 liters/min., contact plate	76076-696
MINI Kit, 200 liters/min., Petri plate	76169-640

MONO




Each kit comes with the following:

TRIO.BAS MONO air sampler, induction battery charger and cable, aspirating head, cover head, ROBUSTUS carrying case, calibration certificate, and fillable IQ, OQ, PQ documents.

VWR Cat. no.

MONO Kit, 100 liters/min., contact plate	76076-704
MONO Kit, 100 liters/min., Petri plate	76169-642
MONO Kit, 200 liters/min., contact plate	76076-706
MONO Kit, 200 liters/min., Petri plate	76169-644

DUO




Each kit comes with the following:

TRIO.BAS DUO air sampler, induction battery charger and cable, aspirating heads, cover heads, ROBUSTUS carrying case, calibration certificate, and fillable IQ, OQ, PQ documents.

VWR Cat. no.

DUO Kit, 100 liters/min., contact plate	76076-708
DUO Kit, 100 liters/min., Petri plate	76169-646
DUO Kit, 200 liters/min., contact plate	76076-710
DUO Kit, 200 liters/min., Petri plate	76169-648

TRIO



Each kit comes with the following:

TRIO.BAS TRIO air sampler, induction battery charger and cable, aspirating heads, cover heads, ROBUSTUS carrying case, calibration certificate, and fillable IQ, OQ, PQ documents.

VWR Cat. no.

TRIO Kit, 100 liters/min., contact plate	76169-650
TRIO Kit, 100 liters/min., Petri plate	76169-654
TRIO Kit, 200 liters/min., contact plate	76169-656
TRIO Kit, 200 liters/min., Petri plate	76169-658

MULTISTATION



Each kit includes:

TRIO.BAS Command Unit Air Sampler, one satellite air sampling chamber with stainless steel aspirating head, cover head, battery charger and cable, cable connection between the satellite and command unit (1 set of 5 meters), LIGHT carrying case (satellite), ROBUSTUS carrying case (command unit), calibration certificate, and fillable IQ, OQ, PQ documents.

VWR Cat. no.

MULTISTATION Kit, 100 liters/min., contact plate **76169-660**

MULTISTATION Kit, 100 liters/min., Petri plate **76169-662**

MULTISTATION Kit, 200 liters/min., contact plate **76169-664**

MULTISTATION Kit, 200 liters/min., Petri plate **76169-666**

Capable of holding up to three satellite units.
Additional satellite kits are listed below.



76169-670

ADDITIONAL SATELLITE KITS with cables

Each kit includes:

One satellite air sampling chamber with stainless steel aspirating head, cover head, cable connection between the satellite and command unit (1 set of 5 meters), LIGHT carrying case.

Contact.....**76169-668**

Petri.....**76169-670**



800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.

ACCESSORIES

DAILY SHIFT HEADS

For use on the MINI, MONO, DUO, TRIO, and Isolator Units. Technopolymer aspirating head. Ready for clean room use. Double bagged and individually packaged. Certificate of sterility assurance provided.

Contact Plate
30/box.....**75804-756**

Petri Plate
30/box.....**75804-758**



75804-758

PORTABLE COMMAND UNIT for data transfer




Portable tablet for remote control of instruments, data download storage and transfer capabilities.
Each.....**76169-672**

CD Software for data transfer
Each.....**75804-738**

VERTICAL HOLDER for hands-free sampling

Allows for easy, hands-free vertical standing of the TRIO.BAS MINI, MONO, DUO, and TRIO air samplers for trolleys or tripods.
Each.....**75804-694**



STAINLESS STEEL VERTICAL HOOK

Sturdy vertical stand for the TRIO.BAS MINI AND MONO.
Each.....**75994-446**
Knob (required)
Each.....**75994-448**



FLOOR TRI-POD

Collapsible stainless steel floor tripod for cleanrooms.
Each.....**75994-454**

