





AC- 5410xx	CUSTOMER DRAWING	Page 1 of 1
Revision 01	Cell Scrapers	 greiner bio-one
	Item-No. 541 0xx	



NOT TO SCALE



M 2:1

Cat.-No.	Description	A	B	C	D
541 070	Cell Scraper	280.0	Ø 7.5	18.0	8.0
541 080	Cell Scraper	400.0	Ø 7.5	18.0	8.0

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date October 24, 2006	Date October 24, 2006	Date October 25, 2006	
Date	Name F. Falter	Name F. Sonntag	Name F. Sonntag	