# 88078 Manganese, AAS standard solution, Specpure®, Mn 1000µg/ml

| | |
|---|---|
| **Product Number:** | 88078 |
| **MDL number:** | MFCD00011111 |
| **Molecular formula:** | Mn in 5% $HNO_3$ |

## Product Specification

| | |
|---|---|
| **Concentration:** | 1,000 µg/ml Manganese |
| **Concentration:** | 1,000µg/mL ±1% |
| **Matrix:** | 5% HNO3 |
| **Starting Material:** | Manganese |
| **Starting Material:** | Mn |

| | |
|---|---|
| **Date of Print:** | January 23, 2022 |
| **Version:** | 1 |

Product Specifications are subject to amendment and may change over time.