PRODUCT DATASHEET

A. PRODUCT DESCRIPTION

Commercial name: Pack of 1000 irradiated microsyringes
Catalog number: 90001780
Catalog dimensions: 320x210x170 mm
Intended purpose: Sterile disposable microsyringes to perform sample spreading using the spiral plater, Eddy Jet. For microbiological use.
Company name: IUL, S.A.

B. SPECIFICATIONS

Working volume: 110µl
Irradiation treatment: Electron-beam

C. PRODUCTS INCLUDED

1. **No. 10101807** 250 irradiated syringes box (4 units). *Figure 1*.





Figure 1. 250 irradiated siringes boxes (4 units).

Each irradiated syringes' box (No. 10101807) includes 25 units of the product **11001941** (10 microsyringes cardboard container). *Figure 2*.

A



B




Figure 2. Closed (A) and opened (B) microsyinges carboard container.

D. IRRADIATED SYRINGES SPECIFICATIONS

Microsyringe. *Figure 3*.

Total dimensions (LxW): 35.7mmx7mm

No. 11001983 - Microsyringe Body: high-density polyethylene copolymer

No. 11003031 - Microsyringe Plunger: polypropylene



Figure 3. Micosyringe body and plunger assembly.

IUL, SA –C/ Ciutat d'Asunción, 4 – 08030 Barcelona Spain **T** +34 93 274 0232 **E** iul@iul-inst.com **W** www.iul-instruments.com