SciLabware Limited
Unit 4, Riverside 2, Campbell Road
Stoke-on-Trent, ST4 4RJ, United Kingdom
Tel +44 (0) 1782 444406
Fax +44 (0) 1782 940436
www.scilabware.com
www.wheaton-uk.com
Registered office (as above)
Registered in England and Wales No. 7723574
BS EN ISO 9001: 2008 registered company

PRODUCT SPECIFICATION

Product Code: 492093

Product Description: 25ml Schellbach Burette, class AS, PTFE key

Brand: PYREX

Type of Glass: Borosilicate glass Type 3.3, which conforms to ISO3585
European and US Pharmacopoeia Type 1 glass

Product Specification: Graduation and inscription in blue enamel
+30s waiting time before final reading can be taken
PTFE key with 1.5mm bore


Grads.	0.10ml
Tol.	±0.05ml
Length	620mm

Product Standard(s): The product is manufactured to comply with ISO385

Pack Quantity: 1 piece per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:



Signed: _____

K. Skupniewicz
Quality Engineer

Date: 13 Feb 2018

Pyrex® is a registered trademark. SciLabware registered user.