Vaerloese, Denmark d. 21-10-2021

Standard Rack for Upright freezers

Product Code: **TE23052**

Side open Rack in stainless steel for upright freezer, for (HxD) 6x3=18 cryoboxes of 32mm

Dimensions (HxDxW):	217x422x139mm
Distance between shelves:	35,50 mm
Maximum box footprint:	137x137 mm
Safety feature:	None
Comfort feature:	Pop-out handle

The Rack is made from high quality Stainless Steel, there is bend and spot welded.

Steel Quality:	EN 1.4307 2R/BA, AISI 304L

The Rack is manufactured in Vaerloese, Denmark.

Best Regards,

TENAK A/S

Erik Abildgaard