BAKERBOND spe, Silica Gel

BAKERBOND spe

Disposable Extraction Columns; Jumbo pack

Material No: 7086-23
Effective from: 2019/11/28

Product Specification

Test	Specification
Average Particle Size, um	about 50
Average Pore Diameter, A	about 60
Loss on drying (at 180°C)	<= 6.0 %
Particle size distribution	35 - 70 um
pH of 10% aquous susp. at 25°C	6 - 8
The following results are based on the silica gel raw material:	
Average Particle Diameter, um (APD)	47 - 60
Specific Surface Area, m2/g	For Information Only
Particle size d90 [um]	<= 105.0
Particle size d10 [um]	>= 33.0

For Laboratory, Research or Manufacturing use.
3mL Solid Phase Extraction Columns
500 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000