2-Propanol

AR

Material No: 6775-45
Effective from: 2013/07/02

Product Specification

Test	Specification
Assay (GC)	>= 99.7 %
Water (KF)	<= 0.1 %
Titrable Acid (meq/g)	<= 0.0005
Titrable Base (meq/g)	<= 0.0001
Residue after Evaporation	<= 0.001 %
Substances insoluble in water	passes test
Aldehydes and ketones (as CO)	<= 0.005 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000