Methyl Alcohol
N.F.

Material No.: 8814-28
:
Revision No.: 0

Product Specification

Meets N.F. Requirements, GMP Manufactured Product

Test	Specification
NF - Assay (CH_3OH) (by GC)	≥ 99.5 %
NF - Identification A	Passes Test
NF - Identification B	Passes Test
NF - Acidity	Passes Test
NF - Alkalinity (as NH_3)	≤ 3 ppm
NF - Nonvolatile Residue	≤ 0.001 %
NF - Readily Carbonizable Substances	Passes Test
NF - Readily Oxidizable Substances	Passes Test
NF - Acetone, Aldehydes (as $(CH_3)_2CO$)	≤ 0.003 %
NF - Water (H_2O)	≤ 0.1 %

GMP Manufactured Product
Bulk Pharmaceutical Chemical
CAUTION: For Manufacturing, processing or repackaging
Appearance (clear, colorless liquid)
No Class 1 or Class 2 Residual Solvents are present in this product.
Only Class 3 Solvents are likely to be present. Residual concentration of all Class 3 Solvents is <0.5%.
Elemental Impurities (USP 232, EP 5.20) - Information on elemental impurities for this product is available on the associated Product Regulatory Data Sheet and elemental impurity profile report.

Packaging Site: Paris Mfg Ctr & DC

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700