REPLIGEN
INSPIRING ADVANCES IN BIOPROCESSING

Repligen Corporation
18617 South Broadwick Street
Rancho Dominguez, California 90220-6435, USA

voice: 310-885-4600
customerserviceUS@repligen.com
www.repligen.com

Date: November 8, 2023

Product Data Sheet

Part Description:	CLOSURE CE 50MM 10/PAK WHITE
Part Number:	142150
Product Type:	Universal Tubing Closure

Specifications, Properties, Recommended Use

Membrane Type:	Standard RC, Biotech RC, Biotech CE	**Function:**	Gentle sealing ridges
Material:	Nylon (polyamide)	**Operating Temperature:**	0 - 90°C
Size:	50mm	**Autoclavable:**	Yes
Color:	White	**Additional:**	n/a
Quantity:	10pc/pk		