GE Healthcare

Automated dissolution sample filtration for increased productivity

Save time and eliminate errors associated with manual sampling by implementing Whatman™ 850-DS 8-channel filter plates in your Agilent™ 850-DS Dissolution Sampling Station*.

Whatman 850-DS 8-channel filter plates developed in conjunction with Agilent offer you the following advantages in your dissolution workflow:

- Whatman 850-DS 8 channel plate allows you to process either 6 or 8 samples at one time economically
- The novel design of the Agilent system and use of the filter plates guaranties 100 percent filtration of samples and reliable alignments of the liquid path reducing the risk of jamming or leaks through the filtration step. The Agilent equipment also allows replacement media to be added to the baths after sample time points to avoid saturation of your dissolution bath.
- Easy to validate the new filter plate format. The filter discs in the plates are made with the same housing and membranes as our 25 mm Puradisc syringe filters.



Ordering information

Description	Product code
850-DS 8-Channel Filter Plate 0.45 µm PTFE 50/PK	7707-3000
850-DS 8-Channel Filter Plate 0.45 µm Nylon 50/PK	7707-3100
850-DS 8-Channel Filter Plate 0.45 µm PES 50/PK	7707-3200
850-DS 8-Channel Filter Plate 0.7 µm GMF 50/PK	7707-3300
850-DS 8-Channel Filter Plate 0.2 µm PTFE 50/PK	7707-3400
850-DS 8-Channel Filter Plate 0.2 µm Nylon 50/PK	7707-3500
850-DS 8-Channel Filter Plate 0.2 µm PES 50/PK	7707-3600
850-DS 8-Channel Filter Plate 0.2 µm PVDF 50/PK	7707-3700
850-DS 8-Channel Filter Plate 0.45 µm PVDF 50/PK	7707-3800

www.gelifesciences.com/Whatman

* Information on the Agilent 850-DS dissolution sampling station can be found at www.agilent.com/lifesciences/850-DS

The image of the Agilent 850-DS Dissolution Sample Station is used with permission. © 2018 Agilent Technologies, Inc.

GE , GE monogram, and Whatman are trademarks of General Electric Company.
Agilent is a trademark of Agilent Technologies, Inc. All other third party trademarks are the property of their respective owner. © 2018 General Electric Company.
GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK.

KA6239181018FL



Agilent 850-DS Dissolution Sampling Station