:\DWGS\CSD-239242


THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF DWK LIFE SCIENCES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES INC. ALL RIGHTS RESERVED.©
ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

TOLERANCE UNLESS SPECIFIED: ±.010

THIRD ANGLE PROJECTION

PRODUCT COMPONENTS

ITEM	QTY	DESCRIPTION



CAP MATERIAL: WHITE POLYPROPYLENE
(F) LINER MATERIAL: .040 WHITE POLYETHYLENE
FOAM WITH .005 NATURAL PTFE-FACED LAMINATED TOGETHER

	TITLE 58-400 CAP, PP WHT, PTFE FCD/FMD PE LNR	DATE DRAWN 06-24-93	PROJ #	–	REVISED TO DWK STDS	07/15/19 EL
		DRAWN BY BILL P. SHIPLEY	APPROX. WGT. 5.5 ±1.5gr	F	2014051410392381	06/05/14 TP
2. APPROXIMATE WEIGHT EXCLUDES LINER.	CAPACITY	CHECKED BY B.c.S.	AUTHORIZED RB	REV	E.C.N. #	DATE
1. DIMENSIONS ARE IN INCHES.	FINISH 58-400	DURAN WHEATON KIMBLE 1501 North 10th Street | Millville, NJ 08332-2038 | USA		DWG. NO: B- CSD-239242		
----NOTES----	SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.					