Technical Data Sheet

Article number	**49102**
Description	**Washing bottle**

Material	LDPE (high pressure polyethylene)
Cap	Screw cap with spray tube
Bottle neck	with DIN thread
Origin	Germany

Art.-Nr.	Volume	Thread
49102018	100 ml	GL 18
49102024	250 ml	GL 25
49102030	500 ml	GL 25
49102036	1000 ml	GL 28

Tel.: +49 9779/808-0
Fax: +49 9779/808-88

Stettener Str. 22-24
97647 Sondheim/Rhön · Germany

E-Mail: info@assistent.eu
Internet: http://www.assistent.eu

Glaswarenfabrik Karl Hecht GmbH & Co KG | Amtsgericht Schweinfurt HRA 8063 | USt.IdNr./VAT-No.: DE132198535
Persönlich haftende Gesellschafterin: Hecht GmbH | AG Schweinfurt HRB1160
Geschäftsführer: Karl-Heinz Hecht, Thomas Diepholtz | Einzelprokura: Tonja Müller