BAKERBOND spe, Octadecyl (C18) Polar Plus/SDB

BAKERBOND spe

Disposable Extraction Columns

Material No: 7704-06
Effective from: 2023/11/08

Product Specification

Test	Specification
Styrene Divinylbenzene (SDB-1):	
Breakthrough Antrazine*	>= 1 L
Breakthrough Desisopropylatrazine*	>= 1 L
Specific Surface Area	>= 950 m2/g
Porosity, Median Pore Diameter	50 - 70 A
Particle Size - 50% Volume	65 - 95
Capacity (mg caffeine / g sorbent)	>= 250
Recovery Atrazine*	>= 97 %
Recovery Desisopropylatrazine*	>= 75 %

For Laboratory, Research or Manufacturing use.
6mL Solid Phase Extraction Columns
Per column: 250 mg Octadecyl (C18) Polar Plus / 100 mg SDB

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000