


UV TUBES

Lamps 6W

Model	Wavelenght (nm)	Lenght (mm)
T-40.L	365	1197
T-40.M	312	1197
T-40.C	254	1197
T-30.L	365	892
T-30.C	254	892
T-20.L	365	586
T-20.C	254	586
T-15.L	365	435
T-15.M	312	435
T-15.C	254	435
T-8.L	365	288
T-8.M	312	288
T-8.C	254	288
T-6.L	365	211
T-6.C	254	211

VILBER
info@vilber.com - www.vilber.com




UV TUBES

Spectral Energy Distribution

365nm






UV-C
Peak wavelength : 365nm
6, 8, 15, 20, 30, 40 watt available
G5 & G13 Cap

312nm






UV-B
Peak wavelength : 312nm
8, 15, 40 watt available
G5 & G13 Cap

254nm





UV-C
Peak wavelength : 254nm
6, 8, 15, 20, 30, 40 watt available
G5 & G13 Cap

Spectral and peak values may differ slightly from the typical values and spectrum above.

VILBER
info@vilber.com - www.vilber.com