Diethyl Ether, Anhydrous

ChromAR HPLC

Stabilized with about 1% ethanol

Material No: 2854-25
Effective from: 2013/06/29

Product Specification

Test	Specification
Assay (GC)	>= 99.5 %
Water (KF)	<= 0.05 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.0005 %
Ultraviolet Absorbance (1.00-cm cell vs water)	
230 nm	<= 0.15
240 nm	<= 0.1
250 nm	<= 0.05
260 nm	<= 0.02
280 nm	<= 0.004
Ethanol	~1 %

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000