







Thread according DIN 168 – GL 45 x 4

Colour	Art.-Nr. (packed)	Art.-Nr. (loose)
Silver	292492807	292492815

Index	Aenderung / changes	Datum / date	Name
D	Knurling has been deleted to make the overcap 292492807 finable	18.12.2019	Schneider
C	Change of design and addjustmnet of dim. after final sampling	05.11.2019	Schneider
B	Change of head plate engaving, add. Art.-Nr.	11.02.2019	Schneider
A	Change of DURAN ® and °C	19.06.2017	Schneider

Weitergabe sowie Vervielfältigung dieses Dokuments, Verwertung und Mitteilung seines Inhalts sind verboten, soweit nicht ausdrücklich gestattet. Zuwiderhandlungen verpflichten zu Schadenersatz. Alle Rechte für den Fall der Patent-, Gebrauchsmuster- oder Geschmacksmustereintragung vorbehalten.
The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

zul. Abweichung
allowed dev.

Oberfläche
surface

	Datum / date	Name
bearb. drawn	07.06.2017	Schneider
geprüft checked	11.02.2019	Rees

Werkstoff / material

Stück quantity

Gewicht weight 0,0045kg

Pos. position

Benennung / title

Aluminium screw closure GL 45
open topped with aperture
Art.-Nr. see table

Maßstab scale 1 : 2

DURAN WHEATON KIMBLE

DURAN Produktions GmbH & Co. KG
www.DWK-LifeSciences.com

Form-Nr. mould number 0GL45

A4

Zeichnungsnummer / drawing number
D366378-GZA-0GL45

Revision D

Blatt/sheet 1 /1

Ers. / repl.f:

Ers. / repl.d:

Schutzvermerk ISO 16016 beachten / Refer to protection notice ISO 16016

INV