Product Data Sheet

Alcoholometers 4 scales

Reference 0500FG000/20-qp



Specifications

Double scale :

- **0-100% vol / gl / Tralles**
  Accuracy : ± 2%
  Div : 1%
- **10-44° Cartier**
  Accuracy : ± 1°
  Div : 0.5°
- Calibration temp. : 20°C
- Reading below meniscus
- ± 290mm
- In plastic box

Alla France
Z.I. du Bompas - 13, rue du Bompas - Chemillé - 49120 Chemillé-en-Anjou - FRANCE
Tél. + 33 (0)2 41 30 55 08 - Fax +33 (0)2 41 30 34 67
www.allafrance.com