Benzyl Alcohol
N.F.

Material No.: 1403-19
Revision No.: 0

Product Specification

Meets N.F. Requirements, GMP Manufactured Product

Test	Specification
NF - Assay ($C_6H_5CH_2OH$)	98.0 - 100.5 %
NF - Acidity	Passes Test
NF - Clarity of Solution	Passes Test
NF - Color of Solution	Passes Test
NF - Identification	Passes Test
NF - Residue on Evaporation	≤ 0.05 %
NF - Peroxide Value	≤ 5
NF - Refractive Index, $n^{20}D$	1.538 - 1.541
NF - Benzaldehyde and Other Related Substances - Benzaldehyde	≤ 0.15 %
NF - Benzaldehyde and Other Related Substances - Cyclohexylmethanol	≤ 0.10 %
NF - Benzaldehyde and Other Related Substances - Total Other Peaks Less Than RT of Benzyl Alcohol	≤ 0.04 %
NF - Benzaldehyde and Other Related Substances - Total Other Peaks Greater Than RT of Benzyl Alcohol	≤ 0.3 %
NF - Benzene	≤ 2 ppm
Solubility	Passes Test

GMP Manufactured Product
Bulk Pharmaceutical Chemical
CAUTION: For Manufacturing, processing or repackaging
The Class 1 residual solvent Benzene is likely to be present below 2 ppm. The actual analysis is reported for each lot.
The
Class 2 solvent Toluene is likely to be present and is process controlled to <30 ppm. No other Class 1, 2 or 3 residual solvents are likely to be present.
This product is not intended for use in preparing parenteral dosage forms.

Packaging Site: Paris Mfg Ctr & DC

James T Ethier

Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700