Technical sheet

Kit Trichrome de Masson, Light Green Variation

Code product: 361350-0000

Composition:

• **Mayer Haemalum**	1 x 100 ml
• **Ponceau Fuchsin**	1 x 100 ml
• **Phosphomolybdic Acid**	1 x 100 ml
• **Light Green**	1 x 100 ml
• **1% Acetic Water**	2 x 100 ml

Applications:

Histo-Cytology:
Staining with kit Trichrome de Masson, Light Green Variation, enables a trichromic staining on histological sections with a specific coloration of the collagen by Light Green.
Mayer Haemalum stains nuclei black-blue or brown.
Ponceau Fuchsin Solution stains cytoplasms red-pink.
Light Green stains collagen green.

Updated: 25/04/2013

RAL Diagnostics - Site Montesquieu - 33651 MARTILLAC Cedex – France - Tel : 33+557 960 404 - Fax : 33+557 960 405 - www.ral-diagnostics.fr