2-Propanol

Baker

Material No: 8119.9025
Effective from: 2017/07/06

Product Specification

Test	Specification
Assay (by GC)	>= 99 %
Acidity (as CH3COOH)	<= 0.002 %
Alkalinity (as NH3)	<= 0.002 %
Boiling Range	<= 1.0 °C
Density (g/ml) at 25°C	0.781 - 0.783
Recorded Boiling Point	82 °C
Residue after Evaporation	<= 0.005 %
Water (H2O)	<= 0.2 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000