Methanol

ChromAR HPLC

Material No: 3041-25
Effective from: 2013/12/05

Product Specification

Test	Specification
Assay (GC)	>= 99.9 %
Water (KF)	<= 0.05 %
Titrable Acid (meq/g)	<= 0.0005
Titrable Base (meq/g)	<= 0.0002
Residue after Evaporation	<= 0.0005 %
Ultraviolet Absorbance (1.00-cm cell vs. water):	
205 nm	<= 1.00
210 nm	<= 0.5
220 nm	<= 0.2
230 nm	<= 0.1
240 nm	<= 0.05
260 nm	<= 0.004
Gradient Elution Test	
235 nm	<= 2 mAU

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000