
DURAN® GLS 80® Laboratory Bottle Wide Mouth Amber

with GLS 80®, USP <660> and USP <671> (Spectral Transmission) compliant





Retrace Code | A 121 °C | USP Standard

With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue quick release closure (PP, integral lip seal) and pouring ring (PP) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: + 140 °C. Alongside easy handling, UV protection up to 500 nm. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology.

Typical applications: storage, transport and safekeeping of light-sensitive substances, easy to use with granulated material, powders and viscous media.

Cat. No.	Capacity (mL)	Thread	d (OD) (mm)	h (mm)	Pack Unit
with screw cap and pouring ring from PP (blue)					
218663653	250	80	95	110	10
1160146	500	80	101	153	10
1160147	1 000	80	101	223	10
1160148	2 000	80	136	253	10
218666959	3 500	80	160	276	1
1160149	5 000	80	182	315	1
1160150	10 000	80	227	390	1
1160151	20 000	80	288	485	1
without screw cap and pouring ring					
218663604	250	80	95	105	10
1178429	500	80	101	148	10
1178430	1 000	80	101	218	10
1178431	2 000	80	136	248	10
218666901	3 500	80	160	271	1
1178432	5 000	80	182	310	1
1178433	10 000	80	227	385	1
1178434	20 000	80	288	480	1