Product Specification Sheet

VWR SLIDE BOXES,PP, FOR 25 SLIDES

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Colour	Pk
631-0738	White	1
631-0739	Blue	1
631-0740	Red	1

Material: **Polypropylene**

W×D×H: **90×140×35 mm**

Other characteristics:

- **Hinged lids are held securely closed by latch lock**
- **Stacking feet maximise stability when placed on top of one another**
- **Slide inventory sheet on the inside of the lid**
- **Friction – fit cover**

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **August 2012**