Nasco SAMPLING
901 Janesville Avenue, Fort Atkinson, WI 53538-0901

Whirl-Pak® Sample Bag
Material Specification Sheet
B01009WA

Size: 3" W x 5" L (7 ½ x 12 ½ cm).

Capacity: 2 oz. (58 ml).

Film: Blend of low-linear and low-density virgin polyethylene, 2.25 mil (.057 mm) thick, FDA approved for food contact.

Closing Tab: Aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak-proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sterilization: Sterilized with ethylene oxide gas, documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality system certified to **ISO9001 & ISO13485.**

Packaging: 500 bags/box in a liner bag twist-tied shut.

Condition of Use: Not to be used at temperatures above 180º F (82º C).