Marienfeld * Am Wöllerspfad 4 * 97922 Lauda-Königshofen * Germany

To whom it may concern

Lauda-Königshofen, 8th January 2021
Ref.: Harry Marienfeld ☎ - 20
E-mail: harry.marienfeld@marienfeld-superior.com

Slide mailer (Art.No. 5600065)

Slide mailers are useful for transport, archiving and staining of microscope slides.



Specifications:

- made of polyethylene (PE-HD)
- break-resistant
- oval shape
- thanks to guide rails inside slides can be accommodated and transported without contacting each other
- with secure closure
- offers space for up to 2 microscope slides (approx. 76 x 26 mm)

Please do not hesitate to contact us should you require any further information. We will be pleased to assist you and look forward to hearing from you.

Since 1922: Our experience is your success.

Best regards,

Harry Marienfeld
President