Honeywell
Specialty Chemicals
Seelze GmbH
P.O. Box 10 02 62
30918 Seelze /Germany
Tel.: 0049 (0)5137-999-0
Fax.: 0049 (0)5137-999-170

SPEZIFIKATION: 02372 N,N-Dimethylacetamide GC-Headspace tested, 99.9%

Print Date: 13.01.2021 / **Revision index:** 03 / **effektiv:** 12.06.2017

suitability for the GC Headspace		complying	
assay (GC)	min.	99.9	%
refractive index (n 20/D)		1.436 - 1.438	
water (Karl Fischer)	max.	0.03	%
absorbance at 268 nm	max.	1.000	
absorbance at 275 nm	max.	0.259	
absorbance at 300 nm	max.	0.070	
absorbance at 350 nm	max.	0.008	
absorbance at 400 nm	max.	0.004	

The above specification shows our present quality of the material.
We reserve the right to change this specification, if necessary.
We warrant the quality of our products according to our General Conditions of Sale.