Chlorobenzene

Baker Analyzed

Material No: 8023.9025ST
Effective from: 2013/05/20

Product Specification

Test	Specification
Assay (by GC)	>= 99 %
Boiling Range	<= 1.0°C
Density (g/ml) at 25°C	1.101-1.104
Recorded Boiling Point	132°C

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000