interscience

BagRack® *Slide*

Storage rack for bags

Ref 221 440



Manufactured under:

MADE IN FRANCE

Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagRack® ***Slide*** enables the clean and methodical storage of 12 bags with samples. Perfect to use with **BagFilter® Pipet** & **Roll**, which allows for an easy pipetting, avoids cross contaminations from the pipet, and allows quick and clean disposal of the bags. The **BagRack®** ***Slide*** allows the user to slide the bag from the side for insertion and also for access for pipetting. This rack is compatible with pre-enrichment of the sample for the detection of pathogens bacteria such as listeria or salmonella.

TECHNICAL SPECIFICATIONS

- Fits 12 bags
- Compatible with:
  - BagFilter® Pipet & Roll for easy pipetting and quick and clean disposal
  - BagFilter® Roll for quick and clean disposal
  - BagClip® and BagFilter® Pipet for easy pipetting
  - BagClip® and standard 400 blender bags
- All stainless steel
- Autoclavable
- Stackable
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

BagRack® *Slide* 400

Ref 221 040

- For 12 x 400 mL bags
- Dimensions (w x d x h): 21 x 35.5 x 22.5 cm
- Box (w x d x h): 23.5 x 37 x 24.5 cm, weight: 1.6 kg

OPTIONS / ACCESSORIES

- **Ref 221 042** **DripTray®:** Stainless steel drip tray for **BagRack®** 400 and **BagRack®** ***Slide*** 400
  Box (w x d x h): 22 x 42 x 4.5 cm, weight: 1 kg

BagRack® slide_FT_1118 Pictures and information are not contractually binding. BagRack® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

www.interscience.com

interscience
30, ch. Bois Arpents - 78860 St Nom - FRANCE
T: +33 (0)1 34 62 62 61
Email: info@interscience.com

interscience USA & CANADA
32 Cummings Park - Woburn, MA 01801 - USA
P: +1 781 937 0007 - F: +1 781 937 0017
Email: sales.usa@interscience.com

interscience CHINA
上海市徐汇区肇嘉浜路798号坤阳大厦1903室-200030
电话: +86 (0)21-64739390 - +86 189 3097 0733
邮 址: sales.china@interscience.com

interscience SOUTH-EAST ASIA
The Centropod, 80 Changi Rd - 05-07 SINGAPORE 419715
T: +65 6909 0825 - M: +65 8118 5178
Email: sales.asia@interscience.com