Seit 1927
über 90 Jahre

Since 1927
more than 90 Years

Amarell
GmbH & Co. KG
Präzisions-Thermometer / Aräometer

Lindenstraße 43
D-97892 Kreuzwertheim

Tel. +49 (0) 9342 9283-0
Fax: +49 (0) 9342 39860

Internet: www.amarell.de
E-Mail: thermometer@amarell.de

Technical Data Sheet

Article number:	M247302-FL/DKD3/CL
Type:	Congealing point thermometer, similar DIN 12785, solid stem, white backed, durable pigment
Measurement range:	+50+150°C
Subdivision:	0,5°C
Filling:	Red special liquid
Length:	300mm
Diameter:	6,5mm
Immersion:	Total
Certificate:	Calibrated, with DKD-Certificate (DAkkS) with 3 test points

September 2021