Product Data sheet

Cytiva Part No.	10531018
Product Description	903 PROTSAVER 5 SPOT 100/PK
Revision:	Rev AD
Release Date:	1-Mar-2021

Physical Parameter:

Specification	Nominal	Min.	Max.
Card Length (fully open) (inches)	4" (±1/16")	3 15/16"	4 1/16"
Card Length (closed) (inches)	2.5" (±1/16")	2 7/16"	2 9/16"
Card Width (inches)	3.5"		
Number of Circles	5		
Perforated circles included? (yes/no)	No		
Shelf Life (years)	5 years		

Technical parameter:

Specification	Nominal	Min.	Max.
Card thickness (mm)	0.108		
Card Weight (grams)	2.95		
Blood Circle Size (mm)	12mm ID		
Blood absorbency (seconds)	8.2	6.8	10.8

Image of Card



