
ThermoFisher SCIENTIFIC

Product Specification

A17682 Tartrazine



Product Number:	A17682
CAS number:	1934-21-0
MDL number:	MFCD00148908
Molecular formula:	$C_{16}H_9N_4Na_3O_9S_2$
Molecular weight:	534.37

Product Specification

Appearance (Color):	Orange
Form:	Powder
Assay (unspecified):	Dye content >85%
Loss on Drying:	≤13.0% (200°C/constant weight)
lambda-max (UV-Vis):	425-430nm (in water)

Date of Print:	June 7, 2022
Version:	3

Product Specifications are subject to amendment and may change over time.