
Unit 14, Graphite Way,
Hadfield, Derbyshire,
SK13 1QH, UK
Tel: +44 (1457) 860 111
Email: enquiries@fluorochem.co.uk
Website: www.fluorochem.co.uk

Product Specification

Product Name: 5-Chlorothiophene-2-carboxylic acid

Catalogue Number	014515
MDL	MFCD00041426
CAS	24065-33-6
Min Assay/Purity	98%
Melting Point	151-155°C
Molecular Formula	C5H3ClO2S
Molecular Weight	162.59



Test	Specification
Appearance	White to Cream to Yellow to Beige to Brown
Physical Property	Crystalline Powder or Solid
Method of Assay	Titration (NaOH)

Fluorochem warrants that, at the time of release, this product should conform to the information contained in this document. For further enquiries please contact customer service. The purchaser must determine the suitability of the product for their own particular use. See website for full terms and conditions.