BAKERBOND spe, Phenyl

BAKERBOND spe

Disposable Extraction Columns

Material No: 7095-03
Effective from: 2017/03/08

Product Specification

Test	Specification
The following results are based on the bulk bonded phase:	
Carbon (C)(%)	For Information Only
Hydrogen (H)(%)	For Information Only
Surface Coverage (ug/m2), Based on %C	185 - 225
The following results are based on the silica gel raw material:	
Particle size distribution < 31 mikrons, wt%	<= 17 %
Particle Size Distribution: APD, microns	47 - 61
Specific Surface Area, m2/g	For Information Only
The following results are based on the solid phase extraction columns:	
Extraction Residue	<= 0.2 %

For Laboratory, Research or Manufacturing use.
3mL Solid Phase Extraction Columns
500 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000