Specification Sheet

SIGMA-ALDRICH

Product Name	Trypan Blue solution, 0.4%, liquid, sterile-filtered, suitable for cell culture
Product Number	T8154
Product Brand	SIGMA
CAS Number	72-57-1
Molecular Weight	960.81

TEST	**SPECIFICATION**
Appearance (Form)	Liquid
Appearance (Colour)	Dark Blue
pH	7.1 - 7.5
Osmolality	260 - 340 mOs/kg
Sterility	Pass
Dye Exclusion (Cell Culture Test)	Pass
Animal Derived Component Free	Conforms