
InLab® Routine Pt1000

pH ● Ref ● °C ● ISM ○

For general laboratory applications; fast and highly chemical resistant

Scale 1:1



Sensor type	combined pH electrode
Order number	51343056

Specifications

pH range	0...14
zero point (pH)	7.00 ± 0.25
temperature range	0...100 °C
temperature probe	Pt1000
shaft material	glass
membrane glass	HA
membrane resistance	< 600 MΩ
reference system	ARGENTHAL™
reference diaphragm	ceramic
reference electrolyte	3 mol/L KCl
bridge electrolyte	-
electrode head	MultiPin™
cable length	-
connector	-
cell constant (cm^{-1})	-
cell type	-
indicating metal	-

Storage

InLab storage solution (30111142)
or 3 mol/L KCl (51350072)

© 12/2015 MTII www.mt.com/electrode-guide

METTLER TOLEDO