Diethyl Ether

AR

Stabilized with BHT

Material No: 6783-25
Effective from: 2013/03/25

# Product Specification

| Test | Specification |
|---|---|
| Assay (GC) | >= 99.5 % |
| Water (KF) | <= 0.2 % |
| Titrable Acid (meq/g) | <= 0.0005 |
| Residue after Evaporation | <= 0.03 % |
| Ethanol | <= 0.1 % |
| Stabiliser (BHT) | ~ 5 ppm |

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000