SIGMA-ALDRICH®


sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Neopentyl glycol diglycidyl ether – technical grade

Product Number:	**338036**
CAS Number:	17557-23-2
MDL:	MFCD00040687
Formula:	C11H20O4
Formula Weight:	216.27 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Miscellaneous Assay	135 - 165
Epoxide Equivalent, g/eq	
Viscosity	10 - 30
mPa.s (25C)	

Specification: PRD.0.ZQ5.10000012925

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.