Type in Product Names, Product Numbers, or CAS Numbers to see suggestions.

Specification Sheet

Product Name	Benzene analytical standard
Product Number	12540
Product Brand	SIAL
CAS Number	71-43-2
Molecular Weight	78.11

TEST	SPECIFICATION
APPEARANCE (COLOR)	Colorless
APPEARANCE (FORM)	Liquid
PURITY (GC AREA %)	≥ 99.9 %
REFRACTIVE INDEX N20/D	1.500 - 1.502
WATER (COULOMETR.)	≤ 0.03 %
PROTON NMR SPECTRUM	CONFORMS TO STRUCTURE
EXPIRY PERIOD	36 MONTHS

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
Reproduction of any materials from the site is strictly forbidden without permission.
Site Use Terms | **Privacy Policy** | **General Terms and Conditions of Sale** | **Copyright Consent** | **Cookies Settings**