Heat Sealer with Cutter

Quickly and efficiently seal polyethylene or polypropylene bags, sheets or flat tubing from 0.01 to 0.8 mm thick



Sliding blade removes excess material 19 mm from the seal

- Can seal a bag up to 12 in (30.5 cm) long
- Heating element measures 31 cm long and is only hot when sealing
- Includes one replacement part kit
- Additional replacement part kits can be ordered separately



Replacement Part Kit:

- 2 extra blades,
- 2 extra heating elements,
- 2 extra PTFE cloths



Timing knob ensures safe sealing

Item No.	Power Supply	L x W X H *in*	L x W X H *cm*
HS28270	220-230 VAC, UK	17.6 x 3.4 x 6.5	44.6 x 8.6 x 16.6
HS28271	220-230 VAC, Cont. Europe/Israel		
HS28272	110-120 VAC, US/Taiwan/Japan		
HS28280	Replacement Part Kit		

Parafilm® M Sealing Film

Protect the integrity of your laboratory and research applications

MADE IN USA

Parafilm® M is a semi-transparent, flexible, thermoplastic film with two-way stretch



- Protect the integrity of research applications with a tight moisture-barrier that holds a strong seal and doesn't release from glassware, despite off-gassing
- Reduce the risk of tampered data with a covering that is odorless, colorless and semi-transparent
- Prevents moisture and volume loss due to material composition and tight seal
- Flexible and self-sealing material allows it to cling around irregular shapes and surfaces
- Stretches more than 200% of original length
- Resistant up to 48 hours against many polar substances such as saline solutions, inorganic acids and alkaline solutions
- Two color options available in 2" rolls (natural and purple) for color coding purposes

Item No.	Color	Film Width *in*	Film Width *cm*	Roll Length *ft*	Roll Length *m*
HS234526A	Natural	2	5	250	76.2
120762	Purple	2	5	250	76.2
HS234526B	Natural	4	10	125	38.1
HS234526C	Natural	4	10	250	76.2
120265	Natural	20	50.8	50	15.2

Parafilm® M is a registered trademark of Amcor Limited