Technical sheet

Kit Gram-Hücker L
slow action

Code product: **362860-0000**

Composition:

• **Crystal Violet Oxalate**	Code: 3614901A240	1 x 240 ml
• **Lugol, PVP-stabilized solution**	Code: 3674002A240	1 x 240 ml
• **Slow Differentiator**	Code: 363030-0240	1 x 240 ml
• **Safranin**	Code: 3615004A240	1 x 240 ml

Applications:

Bacteriology:
Staining with Gram-Hücker L kit enables to differentiate Gram-positive from and Gram-negative bacteria.
Gram-Hücker L kit is perfect for manual and non-simultaneous staining of slides.
Gram-positive bacteria stain violet.
Gram- negative bacteria stain pink.

Parasitology and Mycology:
Staining with Gram-Hücker L kit allows to detect pathogenic fungi as well as microsporidia spores (protozoa).
Pathogenic fungi and microsporidia spores (protozoa) are classified as Gram-positive and stain violet.

Updated: 25/04/2013

RAL Diagnostics - Site Montesquieu - 33651 MARTILLAC Cedex – France - Tel : 33+557 960 404 - Fax : 33+557 960 405 - www.ral-diagnostics.fr