
| CATALOG NO. | NOMINAL CAPACITY, (GAL.) | | MATERIALS | |
|---|---|---|---|---|
| | | | BOTTLE | CLOSURE |
| 002120 | 1/2 | 1 | HDPE | PP |
| 002121 | 1/2 | 1 | PP | PP |

| SIZE CODE.....: | -0005 | | -0010 | |
|---|---|---|---|---|
| SIZE (NOM VOL L/GAL) | 2/ 1/2 | | 4/ 1 | |
| NOM. BRIM CAP. (mL): | 2200 | | 4300 | |
| NECK FINISH (mm)...: | 100-415 | | 100-415 | |
| THREAD (T.P.I.)....: | 5 | | 5 | |
| NOM. WEIGHT W/CLOSURE (gms)....: | 200 | | 300 | |
| | IN. | MM | IN. | MM |
| A (NECK I.D.).....: | 3.51 ±.08 | 89.2 ±2.0 | 3.51 ±.08 | 89.2 ±2.0 |
| B (HT. W/CLOSURE)..: | 9.60 ±.12 | 243.8 ±3.0 | 11.32 ±.12 | 287.5 ±3.0 |
| C (HT. W/O CLOSURE): | 9.26 ±.12 | 235.2 ±3.0 | 10.98 ±.12 | 278.9 ±3.0 |
| D (DIA.).....: | 4.73 ±.09 | 120.1 ±2.3 | 6.02 ±.10 | 152.9 ±2.5 |
| W.T.(NOM. WALL THK): | 0.080 | 2.03 | 0.080 | 2.03 |




NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.


CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.



Therno Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA


| MADE MG | APP'D PVC 2/25/92 | PSS 00-212X-XXXX-1 | REV. 0292 |
|---|---|---|---|
| CK'D PVC | | | |

DESCRIPTION:

# LARGE WIDE-MOUTH BOTTLES