Specification Sheet

Product Name	CDI ≥97.0% (T)
Product Number	21860
Product Brand	ALDRICH
CAS Number	530-62-1
Molecular Weight	162.15

TEST	SPECIFICATION
APPEARANCE (COLOR)	White to Faint Yellow
APPEARANCE (FORM)	Powder or Crystals
TITRATION	97.0 - 103.0 %
TITRATION (METHOD)	ACIDIMETRIC
MELTING POINT	115 - 122 C
PROTON NMR SPECTRUM	CONFORMS TO STRUCTURE
RECOMMENDED RETEST PERIOD	36 MONTHS

Feedback