Technical Sheet

FIX-RAL 555

Product code: 362870

Specifications:

- **Appearance:** Clear violet Solution
- **Absorption Maximum in Methanol:** 586 nm ≤ λ ≤ 588 nm

Applications:

Haematology: used in Kit RAL 555, fast variation of May-Grünwald Giemsa technique, for differential stainings of medullary and blood smears.

Parasitology and Mycology : used in Kit RAL 555 for the detection of parasites and fungi on blood smears.

Histo-Cytology: used in Kit RAL 555 for polychromic quick stainings of histological sections and the detection of Helicobacter pylori on histological sections.

Updated: 07/02/2022