HIRSCHMANN®

Data sheet

Article: Automatic burettes acc. to Dr. Schilling 5 ml
Code no.: 360 01 53

Material:
DURAN (Borosilicat glass 3.3 according ISO 3585)

Chemical properties:
- Hydrolytic resistance at 98°C: Class 1 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 1 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

AmberTube:
Natural rubber (NR)

Discharge tip clear:
Polyvinylchlorid (PVC)

Bottle:
Low-density polyethylene (LDPE)

Nominal volume:
5 ml

Subdivision:
0,05 ml

Tolerance:
± 0,03 ml

Standard:
ISO 385

Hirschmann Laborgeräte
GmbH & Co.KG

Roland Benz
Quality Manager

Eberstadt, 25.07.2022